

Autoclave Vertical



MSIW- 034

SPECIFICATIONS & FEATURES:-

- **Autoclave (Vertical, Automatic) with Digital Temperature Indicator** – suitable for sterilizing and discarding.
- **Material:** Complete Stainless-Steel body (SS 304 quality).
- **Construction:** Double-walled design with Boiler, Outer Shell, and Lid – all stainless steel.
- **Temperature Monitoring:** Digital temperature indicator.
- **Safety & Operation Features:**
 - Radial locking system with paddle lifting device.
 - Automatic water level cut-off device.
 - Automatic pressure control switch.
 - Neoprene rubber gasket (joint-less).
 - ISI-marked immersion heating element.
- **Performance:**
 - Pressure ranges up to 15 psi.
 - Capacity: **35 Liters.**
 - Size: **250 mm × 400 mm.**
- **Electrical Requirements:**
 - Operates on **220V AC, Single Phase.**
 - Rating: **1.5 kW.**
- **Supplied Complete:** With cord and plug.

M. S. Scientific Instrument Works

Country of Origin : INDIA

E-112, KH. NO. 77/9, BLK-E, RAJEEV NAGAR EXT.,
NEAR RAHUL PUBLIC SCHOOL, NEW DELHI- 110080 (INDIA)

