

HOT PLATE



MSIW- 016

Features & MOC	For continuous heating up to 250 deg C. the body is made of specially designed perforated P.C.R.C. Sheet painted with attractive stove enamel with insulated elements with Kanthal wire.
Size	8'' Dia
Power Supply	To work on 220/230 Volt s.

M. S. Scientific Instrument Works

Country of Origin : INDIA

E-112, KH. NO. 77/9, BLK-E, RAJEEV NAGAR EXT.,
NEAR RAHUL PUBLIC SCHOOL, NEW DELHI- 110080 (INDIA)

