

Rectangular Water Bath



MSIW- 017

SPECIFICATIONS & FEATURES :-

Construction:

- **Double-walled design:**

Inner Chamber: Stainless Steel (SS 304 grade).

Outer Body: Stainless Steel (SS 304) with attractive stove-enamel finish.

- **Fitted with concentric rings** for uniform heating.
- **On/Off switch** provided for easy operation.

Temperature Control:

- Controlled by **hydraulic thermostat**.
- **Temperature Range:** 5°C above ambient to 100°C.
- **Sensitivity/Accuracy:** $\pm 0.5^{\circ}\text{C}$.

Power Requirements:

- Operates on **220/230 V AC**.

Capacity/Size-300 X 250 X 100mm (6 of 75mm Dia)

M. S. Scientific Instrument Works

Country of Origin : INDIA

E-112, KH. NO. 77/9, BLK-E, RAJEEV NAGAR EXT.,
NEAR RAHUL PUBLIC SCHOOL, NEW DELHI- 110080 (INDIA)

