

Soxhlet Extraction Apparatus



MSIW-054

SPECIFICATION & FEATURES :

Features & Construction:

- Suitable for solvent extraction in laboratory applications.
- Provided with **three hot plates** (individual operation).
- **Temperature Control:** Each hot plate fitted with individual **Sunvic energy regulator** and indicator light.
- **Maximum Surface Temperature:** Up to **350°C**.
- **Body:** Mild Steel construction, durable and powder-coated.
- **Support:** Two horizontal and vertical support rods with gutkas (clamps) for flask holding.

Flask Capacity:

- Compatible with flasks from **50 ml to 500 ml**.

Power Requirements:

- Operates on **220/230 V AC supply**.

M. S. Scientific Instrument Works

Country of Origin : INDIA

E-112, KH. NO. 77/9, BLK-E, RAJEEV NAGAR EXT.,
NEAR RAHUL PUBLIC SCHOOL, NEW DELHI- 110080 (INDIA)

